[Gebrauchsanleitung](#) | [Operating manual](#) | [Mode d'emploi](#) |
[Instrucciones de manejo](#) | [Istruzione](#) | [Instruções de utilização](#) | [操作手册](#) |
[Руководство по эксплуатации](#) | [사용지침](#) | [Bruksanvisning](#) | [使用説明書](#) |
[Használati utasítás](#) | [Návod k použití](#) | [Gebruiksaanwijzing](#) |
[Instrukcja użytkowania](#) | [Kullanım Talimatları](#)

Transferpette® S -8/-12

Pipettierhelfer | Pipette controllers

Impressum

BRAND GMBH + CO KG
Otto-Schott-Str. 25
97877 Wertheim (Germany)

T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

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Indien:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791
info@brand.co.in
www.brand.co.in

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Table of contents

1 Introduction	31
1.1 Scope of delivery	31
1.2 Use the instructions for use.....	31
2 Safety regulations.....	32
2.1 Safety regulations	32
2.2 Purpose	32
2.3 Limitations of use.....	32
2.4 Application restrictions	33
2.5 Operating exclusions	33
3 Functional and operational components	34
4 Pipetting	35
5 Checking the volume	38
6 Accuracy table	39
7 Adjustment – Easy Calibration.....	40
8 Disinfection/autoclaving.....	42
8.1 Autoclaving.....	42
8.2 UV sterilization	42
9 Maintenance.....	43
9.1 Separate the pipetting unit from the handle	43
9.2 Disassembly/cleaning.....	43
10 Troubleshooting	49
11 Product markings	50
12 Order Information.....	51
12.1 Ordering Information/Accessories..	51
12.2 Spares	52
12.3 Additional accessories	52
13 Repairs.....	53
13.1 Sending for repair	53
14 Calibration service	55
15 Warranty	56
16 Disposal	57

1 Introduction

1.1 Scope of delivery

BOECO - micropipette, with quality certificate, 2 x TipBox filled with BRAND pipette tips, shelf holder, reagent reservoir, mounting wrench, silicone grease and 1 set of shaft seals

1.2 Use the instructions for use

- Carefully read the operating manual before using the device for the first time.
- Keep the user manual easily accessible. It is part of the device.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find up-to-date versions of the operating manual on our website: www.brand.de.

1.2.1 Signal words and their meaning

Signal word	Meaning
⚠ WARNING or ⚠ WARNING! ...	WARNING indicates a dangerous situation that, if not avoided, could result in death or serious injury.
⚠ CAUTION or ⚠ CAUTION! ...	CAUTION indicates a hazardous situation that, if not avoided, could result in moderate or minor injury.
NOTICE or NOTICE ...	NOTE is used to address actions that are not related to physical injury. Example: Possible property damage.

1.2.2 Presentation of descriptions of actions

Format	Meaning
1. Task	Indicates a task.
a., b., c.	Indicates the individual steps of a task.
>	Indicates a prerequisite for a task.
⇒	Indicates a result of a completed task.

2 Safety regulations

2.1 Safety regulations

Please read carefully!

The laboratory device BOECO - micropipette can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

- Every user must read and observe the user manual enclosed with this laboratory device before and during use. The device may only be used by trained and qualified personnel.
- Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.
- When working with infectious or hazardous samples/media (e.g. hazardous substances), the general safety rules in the laboratory must be observed and regulations for handling the samples/media must be followed. The information provided by the media manufacturer (e.g. safety data sheets) must be observed.
- Use the laboratory device only for pipetting or dispensing media within the defined limitations and restrictions of use. Comply with the operating exclusions.
- If working with inflammable media, make sure to take the necessary precautions to avoid the buildup of static electricity (e.g., do not dispense into plastic vessels and do not wipe instruments with a dry cloth). Do not use the device in potentially explosive atmospheres. If in doubt, contact the manufacturer or supplier.
- Always check that the instrument is in proper working condition before use. If device malfunctions are indicated (e.g. sluggish pistons, leaks or power supply faults), stop working with the device immediately and refer to the "Troubleshooting" section in the user manual. Contact the manufacturer, if necessary.
- Always perform work in a manner that does not endanger users or other people. Avoid splashes. Only use suitable vessels. Never use unnecessary force when operating, cleaning or servicing the laboratory device.
- If the laboratory device is supplied with power through a power adapter, batteries or accumulators, the components and connections on the device must be regularly checked for proper working condition. Do not operate the laboratory device and its accessories in an unprotected, humid or wet environment.
- Do not make any technical modifications. Use only original replacement parts of the manufacturer; even power adapters or accumulators of identical size and specification from other manufacturers are not permitted. Do not disassemble the laboratory device and its accessories (e.g., power adapters, cables, racks, accumulators or batteries) any further than described in the user manual!
- Only autoclave the laboratory device when this is permissible according to the user manual.

2.2 Purpose

This is an air displacement pipette for pipetting liquids of low to medium viscosity.

2.3 Limitations of use

The pipette is used for dispensing liquids within the following limits:

- Operating temperature from +15 °C to +40 °C (59 °F to 104 °F). Additional temperatures upon request
- vapor pressure up to 500 mbar

- Viscosity: 260 mPa s

With viscous liquids, the dispensing speed may need to be adjusted.

2.4 Application restrictions

Viscous and wetting liquids may compromise volumetric accuracy. Volumetric accuracy may also be affected when pipetting liquids whose temperature deviates from the ambient temperature by more than $\pm 1\text{ }^{\circ}\text{C}/\pm 1.8\text{ }^{\circ}\text{F}$.


2.5 Operating exclusions

The user must verify the suitability of the instrument for the intended purpose, as aggressive liquids and their vapors can damage the instrument (corrosion!). The instrument cannot be used for the following liquids:

- Liquids with very high steam pressure
- Liquids that corrode the following materials
 - Fluoroelastomer rubber (FKM)
 - Polycarbonate (PC)
 - Polyphenylsulphide (PVDF) (with 50 μl , 100 μl , 200 μl and 300 μl instruments)
 - Polypropylene (PP)
 - Polyvinylidene fluoride (PVDF)
 - Silicone

Additional information on the chemical resistance of plastics can be found at www.brand.de.

3 Functional and operational components

- | | |
|------------------------------------|-------------------------------|
| 1 Pipetting button | 2 Volume-setting wheel |
| 3 Easy Calibration function | 4 Finger rest |
| 5 Pipetting unit | 6 Tip cone |
| 7 Handle | 8 Volume display |
| 9 Volume-change protection | 10 Tip ejection key |

Label window

The instrument can be individually labeled on the finger rest:

- a. Remove the label window on the finger rest.
- b. Mark the labeling film.
- c. Reinsert the labeling film with window.


4 Pipetting

1. Inserting a tip

NOTICE


Perfect analysis results can only be achieved by using quality tips. We recommend pipette tips from For additional information, refer to the Accuracy table, p. 39.

Pipette tips are disposable products!

- a. Insert tip vertically:
Use the correct tips, in accordance with the volume range or color code!
Make sure that the tips are firmly in place and leak tight.

2. Setting the volume

- a. Slide the volume-change protection upward (UNLOCK).
- b. Turn the volume-setting wheel to select the desired volume. In doing so, turn the adjustment wheel steadily, avoiding abrupt turning motions.
- c. Slide the volume-change protection downward (LOCK). The volume-setting wheel becomes noticeably more difficult to turn, but movement is not completely blocked!



3. Align pipetting unit

The pipetting unit can turn freely in both directions.

4. Aspirating a sample

NOTICE

The ISO 8655 standard requires that pipette tips are pre-wetted once before the actual pipetting procedure.

- c.** Allow the pipetting button to steadily move back to its original position.


Leave the tip immersed in the liquid for a few seconds, so that the set volume is aspirated completely. This is particularly important for viscous media and for pipettes with a large volume.

Volume range	Immersion depth	Wait time
0.5 – 100 µl	2 – 3 mm	1 s
100 – 300 µl	2 – 4 mm	1 s
> 1000 µl	3 – 6 mm	3 s

NOTICE

Do not lay the instrument down when the tip is filled; this can cause the medium to flow into the instrument and contaminate it! The instrument should always be stored in the provided shelf mount or table stand and kept in an upright position, without any tip inserted.

5. Dispensing a sample

- a.** Place the pipette tip against the vessel wall. Hold the pipette at an angle of 30-45° to the vessel wall.
- b.** Press the pipetting knob at a steady speed until it stops and hold it down. To improve accuracy, comply with the corresponding wait time for serums, highly-viscous or low-density media.


- c. Completely empty the tip by over-stroking: Press the pipetting button as far as it will go.
- d. While doing this, wipe the pipette tip against the vessel wall.
- e. Remove the pipette tip from the vessel wall and allow the pipetting button to move back to its original position.

6. Ejecting a tip

NOTICE

The instrument should always be stored in the provided shelf mount or table stand and kept in an upright position, without any tip inserted.

- a. Hold the pipette shaft over a suitable disposal container and press the tip eject button down as far as it will go.

5 Checking the volume

We recommend testing the instrument every 3 to 12 months depending on the level of use. However, the testing cycle can be adapted to meet individual requirements. The complete testing procedure (SOP) can be downloaded at www.brand.de.

You can download the detailed test instructions (SOP) at www.brand.de. For GLP- and ISO-compliant evaluations and documentation, we recommend the EASYCAL™ calibration software from BRAND. A demo version can be downloaded from <https://shop.brand.de/>.

Gravimetric volume testing of the pipette is carried out according to the following steps and complies with DIN EN ISO 8655:2022.

- a. Set the maximum specified instrument volume (for procedure, see Pipetting, p. 35).
- b. Condition the pipette before testing by aspirating and dispensing the test liquid (distilled water) with a pipette tip five times.
- c. Aspirate the test liquid and pipette into the weighing vessel.

NOTICE

Each individual channel must be inspected separately.

- d. Weigh the pipetted amount with an analysis scale. (refer to the operating manual of the balance manufacturer.)
- e. Calculate the pipetted volume. In doing so, take into account the temperature of the test liquid.
- f. At least 10 pipetting series and weighings in three volume ranges (100%, 50%, 10%) are recommended. Two tips must be used for each volume range to be tested.

Calculation (for nominal volume)

x_i = weighing results

n = number of weighings

V_0 = nominal volume

Z = Correction factor (e.g. 1.0029 µl/mg at 20°C, 1013 hPa)

Mean:

Mean volume:

Accuracy*:

$$\bar{x} = \frac{\sum x_i}{n}$$

$$\bar{V} = \bar{x} * Z$$

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Standard deviation*:

Coefficient of variation*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

$$CV\% = \frac{100 s}{\bar{V}}$$

*) Accuracy and coefficient of variation are calculated according to the formulas of statistical quality control.

6 Accuracy table

Volume range [μl]	Partial volume [μl]	A* ≤ ±	CV* ≤%	Sub steps [μl]	Recommended tip type [μl]
0.5–10	10 5 1	1.6 2 8	1.0 2 6	0.01	0.5...20
5...50	50 25 5	0.8 1.4 6	0.4 0.8 3	0.05	2...200
10...100	100 50 10	0.8 1.4 4	0.3 0.6 2	0.1	2...200
20...200	200 100 20	0.8 1.4 4	0.3 0.6 1.5	0.2	2...200
30...300	300 150 30	0.8 1.4 4	0.3 0.6 1.5	0.5	5...300


*A = Accuracy, CV = Coefficient of Variation


Final test values based on the nominal volume (= max. volume) printed on the device and the specified partial volumes at the same temperature (20 °C/68 °F) of the device, surroundings and distilled water, in accordance with DIN EN ISO 8655.

7 Adjustment – Easy Calibration


The instrument is permanently calibrated for aqueous solutions. If it is determined that the pipette is operating inaccurately or to adjust the instrument to work with solutions of varying density and viscosity or with specially-shaped pipette tips, it can be calibrated using the Easy Calibration Technique.


- a. Perform a volume check and determine the actual value; see Checking the volume.
- b. Remove label window (1) and labeling film (2): Gently move the clamp and lift it off.
- c. Using a paper clip or an unused pipette tip, remove the protective film (3) (the protective film can be discarded).

- d. Slide the red adjustment slider back completely, lift the volume-setting wheel (decoupling) and release the adjustment slider.

- e. Set the adjustment value:
Transferpette® S, adjustable: with the volume-setting wheel in the UNLOCK position, set to the previously determined actual value.
Transferpette® S, fixed-volume: set the volume by rotating in the +/- direction.
 A volume check is recommended after every adjustment.

- f. Slide the adjustment slider completely back again, push the volume-setting wheel downward and release the adjustment slider. Re-attach the labeling film and re-assemble the label window.

NOTICE

The change to factory settings is indicated by the red adjustment slider now visible in the label window.

8 Disinfection/autoclaving

8.1 Autoclaving

The Pipette is completely autoclavable at 121 °C (250 °F), 2 bar and a holding time of at least 15 minutes, in accordance with DIN EN 285.

- a. Eject the pipette tip.
- b. Autoclave the complete pipette without any further disassembling.
- c. Allow the pipette to completely cool and dry.

NOTICE

The effectiveness of autoclaving must be verified by the user. Maximum safety is achieved through vacuum sterilization. We recommend the use of sterilization bags.

NOTICE

Prior to autoclaving, the volume-setting wheel must be set on an available numbered value (e.g., 11.25 or 11.26, but not between), with the volume-change protection set to fully unlocked (UNLOCK).

If the pipette is autoclaved frequently, the piston and seal should be greased with the supplied silicone grease in order to ensure proper movement. After autoclaving, tighten the connection between the hand grip and the pipette shaft if necessary.

8.2 UV sterilization


The device is resistant to normal exposure to a UV disinfection lamp. The effects of the UV exposure may cause some color change.

9 Maintenance

- a. Check the pipette tip cone for damage.
- b. Inspect the piston and seal for contamination.
- c. Check the device for leaks.

We recommend using the BRAND leak detector, the BRAND PLT unit. As an alternative to this, aspirate a sample and hold the device vertically for approx. 10 s. If a drop forms at the pipette tip, follow the malfunction remedy, see Troubleshooting, p. 49.

9.1 Separate the pipetting unit from the handle

- a. Eject the pipette tips.
- b. To separate, pull the pipetting unit downward as far as possible. **Only then**, turn it clockwise. After one revolution, the unit no longer needs to be pulled downward while turning.

NOTICE


- > When assembling, the pipetting unit must be screwed counter-clockwise onto the handle so that it clicks into place.
- > When assembling, do not pull the pipetting unit downwards.
- > Improper handling can lead to damage.

9.2 Disassembly/cleaning

The three main components of the pipetting unit can be easily separated and disassembled for maintenance, cleaning or replacing parts.

Changing the O-rings on the individual shafts is described in detail in the instructions included with the replacement part.

Main components of the pipetting unit

A

Piston unit with piston support bar [I] and inserted pistons, which can be individually unscrewed for cleaning or replacement.

B

Shaft unit with shaft support bar [II] and attached central guide axis (Z) as well as the shafts and seals, which can be individually unscrewed for cleaning and replacement.

C

Pipette housing, which is connected to the pipette housing cover [III] of the piston unit using two turn-lock closures.

Maintenance

In order to ensure proper functioning, the Transferpette® S -8/-12 should be serviced at regular intervals and cleaned as necessary.

What must be checked?


- a. Check the pipette shafts, pistons and seals for damage and contamination.
- b. Check the device for leaks. We recommend using BRAND's leak detector, the BRAND PLT unit. As an alternative to this, aspirate a sample and hold the device vertically for approx. 10 s. If drops form on the pipette tips, refer to "Troubleshooting" Troubleshooting, p. 49.

NOTICE


Cleaning instructions

- > Clean the individual shafts and pistons (only these parts) with a soap solution or isopropanol, and then rinse with distilled water.
- > Allow the parts to completely dry and cool. Liquid residues in the shafts result in deviations in accuracy.
- > Re-grease pistons with a very thin coat of the supplied grease. For the central guide axis (Z), use only the prescribed fluorine static grease!


Removing shafts and seals for cleaning or replacement


- a. Separate the pipetting unit from the hand grip.
- b. Rotate both closures of the pipette housing cover by 90° (e.g. using a coin) and pull off the pipette housing.


- c. Place the installation tool (M) onto a single shaft and unscrew the shaft.

- d. Push the piston unit all the way down. After removing the shaft, the seal is located either in the shaft or on the piston. Remove and inspect the seal and clean or change as necessary. If required, re-grease the piston with the supplied silicone grease. (To clean the Transferpette® S -8/-12 30 - 300 µl, remove the additional pressure ring from the piston.)


- e. Slide the seal onto the piston with the flat side upwards. (On the Transferpette® S -8/-12 30 - 300 µl, first re-install the additional pressure ring!) Tighten the cleaned or new shaft using the installation tool.


- f. Reassemble the pipetting unit. Screw the pipetting unit counterclockwise into the handle. Gently press the pipetting unit into the handle while screwing it in. Then check the instrument for leaks, see Maintenance, p. 43.

Removing pistons for cleaning or replacement


- a. Separate the pipetting unit from the hand grip.
b. Rotate both closures of the pipette housing cover by 90° (e.g. using a coin) and pull off the pipette housing.


- c. Remove both outer Phillips-head screws on the shaft unit.

NOTICE


The central guide axis (Z) must not be removed!


- d.** Pull apart and separate the piston and shaft unit. Remove stroke springs.


- e.** Place the installation tool onto the piston nut and unscrew it.

- f.** Remove the piston nut and pull off the piston with piston spring.

- g.** Insert piston spring (a) and cleaned or new piston. Retighten the piston nut using the installation tool. Insert stroke springs (b).


- h.** Reassemble the pipetting unit. Screw the pipetting unit counterclockwise into the handle. Gently press the pipetting unit into the handle while screwing it in. Then check the instrument for leaks, see Maintenance, p. 43.

10 Troubleshooting

Problem	Possible cause	Corrective action
Tip dripping (device leaking)	Unsuitable tip	Only use high-quality tips
	Tip not seated tightly	Firmly press tip on
The instrument does not aspirate or aspirates too little; the dispensed volume is too low	Seal contaminated	Clean seal
	Seal or cone is damaged	Replace seal or shaft
	Piston is contaminated or damaged	Clean or replace piston
Aspiration is very slow	Shaft is clogged	Clean shaft
Dispensed volume too large	Pipetting button pressed too far (to the over-stroke point) before aspirating	Ensure proper handling.
Piston sluggish	Piston is contaminated or not greased	Clean piston and apply grease

11 Product markings

Symbol or number	Meaning
	Read the user manual.
XXZXXXXX	Serial number
	The instrument is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation. Character sequence DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was affixed.
	Autoclavable up to the temperature shown
Data matrix code or quick response code	The codes link to the BRAND MyProduct website.
www.brand.de/ip	Hyperlink to BRAND patent site

12 Order Information

12.1 Ordering Information/Accessories

Transferpette® S -8

Volume	Description	Order No.
0.5–10 µl	M8-10	705900
5–50 µl	M8-50	705906
10–100 µl	M8-100	705908
20–200 µl	M8-200	705910
30–300 µl	M8-300	705912

Transferpette® S -12


Volume	Description	Order No.
0.5–10 µl	M12-10	705920
5–50 µl	M12-50	705926
10–100 µl	M12-100	705928
20–200 µl	M12-200	705930
30–300 µl	M12-300	705932

Accessories

	Description	Order No.
	Table stand for 6 single-channel or multi-channel pipettes.	704807
	Table stand for 1 single-channel or multi-channel pipette.	703440
	Wall mount	704812
	Shelf/rack mount	704811

12.2 Spares

12.2.1 Volumes up to 300 µl

Volume [µl]	Piston	Shaft ^{1,2}	Seal	O-ring
0.5-10 µl	705659	705677	703340	703380
5-50 µl	705666	705634	703343	705618
10-100 µl	705662	705635	703344	705618
20-200 µl	705663	705636	703345	705618
30-300 µl	705664	705638	703346	705618

¹Until March 23, other order no.; see shop.brand.de

²incl. seal, O-ring, and installation tool. BOECO - micropipette 30–300 µl also with pressure ring.

NOTICE The appearance and dimensions of the spare parts correspond to the respective nominal volume.

12.3 Additional accessories

Description	Order No.
Label window, PU 1 pc.	704750
Labeling film, PU 5 pcs.	704751
Silicone grease	703677
Fluorine static grease	703678
Reagent reservoir, PP. Content 60 ml. Autoclavable at 121 °C (2 bar), with lid. PU 10 pcs.	703459
Reagent reservoir, PP. Content 60 ml. Autoclavable at 121 °C (2 bar), sterile, without lid. Packaged individually. PU 100 pcs.	703411
Reagent reservoir, PP. Content 60 ml. Autoclavable at 121 °C (2 bar), sterile, without lid. 5 pcs./bag. PU 200 pcs.	703409
PLT unit (pipette leak detector)	703970

13 Repairs

13.1 Sending for repair

NOTICE

Transporting hazardous materials without approval is prohibited by law.

Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Within USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

Germany:

BRAND GMBH + CO KG
Otto-Schott-Straße 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T +1-860-767 2562
F +1 - 860 - 767 2563
info@brandtech.com
www.brandtech.com

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791
info@brand.co.in

China:

BRAND (Shanghai) Trading Co., Ltd.
Rm 201-202, North Tower,
No. 199 Kaibin Rd, Xuhui District, Shanghai
Shanghai 200030 (P.R. China)
T +86 21 6422 2318
F +86 21 6422 2268
info@brand.com.cn
www.brand.cn.com

14 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the device. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive Liquids.

The complete SOP for testing can be downloaded from www.brand.de or www.brandtech.com.

BRAND also offers you the option of having your devices calibrated through our factory calibration service or through our accredited calibration laboratory. Just send us the devices to be calibrated, indicating the type of calibration you would like. You will get your devices back in a few days. A detailed calibration report (factory calibration) or an accredited calibration certificate in accordance with DIN EN ISO/IEC 17025 is enclosed with each device. More information can be obtained from your retailer or directly from BRAND. The order document is available for download at www.brand.de (Service & Support).

For customers outside Germany

If you would like to use our calibration service, please contact one of our service partners in your region. Our service partners can forward your devices to BRAND for factory calibration, if required.

15 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the instrument or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

16 Disposal

Before disposal, observe the relevant national disposal regulations and ensure that the product is disposed of properly.

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