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## Notes prior to the use of the manometers


The manometers used are vacuum meters or overpressure manometers with glycerin filling. The overpressure manometer indicates the overpressure relative to the atmospheric pressure at the place of installation. The vacuum meter indicates the absolute pressure.

- ➔ Install and tighten the manometers first (see instructions for use of the pump) or install the pump with manometer in the room of operation.
- ➔ Prior to use the reference chambers of some of the manometers (depending on the model of manometer) have to be ventilated.
- ☞ Else, not ventilating the manometer may lead to a systematic measuring error.

There are different types of manometers. To ventilate the reference chamber check which types of manometers are installed at the pump and proceed according to the given instructions:


pressure compensation valve

### A: Manometer with pressure compensation valve

Ventilate the manometer using the pressure compensation valve. Close the valve after ventilation.


- ☞ Repeat if necessary until dial shows zero with measuring connection ventilated to atmosphere.

pin

### B: Manometer with pin

Ventilate the manometer by lifting the pin. Do not remove the pin completely!

hole

### C: Manometer with bore at the rear side

Pierce the hole at the rear side of the manometer. Some outflow of liquid is normal.

- ☞ Do not tilt pump during transport.

### D: Manometer without bore, pin or pressure compensation valve

Do **not** ventilate the reference chamber of this manometer!

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