




PCR-Platte 96-well, rigid frame PC/PP, frosted, wells transparent 781541






L	127.76 mm
W	85.48 mm
H	16.1 mm
H _w	15.1 mm
X1	14.38 mm
X _n	113.38 mm
Y1	11.24 mm
Y _n	74.24 mm
P	9.00 mm
D	5.46 mm
C	12-well
R	8-well

Well Number	96
Type	full skirt, low profile
Nominal Volume	150 µl
Material	PC/PP
Cut corner	H1
Quality	BIO-CERT® PCR QUALITY



Subject to technical modification without notice. Errors excepted.

PCR-Platte 96-well, rigid frame PC/PP, golden, wells transparent 781542






L	127.76 mm
W	85.48 mm
H	16.1 mm
H _w	15.1 mm
X1	14.38 mm
X _n	113.38 mm
Y1	11.24 mm
Y _n	74.24 mm
P	9.00 mm
D	5.46 mm
C	12-well
R	8-well

Well Number	96
Type	full skirt, low profile
Nominal Volume	150 µl
Material	PC/PP
Cut corner	H1
Quality	BIO-CERT® PCR QUALITY



Subject to technical modification without notice. Errors excepted.

PCR-Platte 96-well, rigid frame PC/PP, red, wells transparent 781543






L	127.76 mm
W	85.48 mm
H	16.1 mm
H _w	15.1 mm
X1	14.38 mm
Xn	113.38 mm
Y1	11.24 mm
Yn	74.24 mm
P	9.00 mm
D	5.46 mm
C	12-well
R	8-well

Well Number	96
Type	full skirt, low profile
Nominal Volume	150 µl
Material	PC/PP
Cut corner	H1
Quality	BIO-CERT® PCR QUALITY



Subject to technical modification without notice. Errors excepted.

PCR-Platte 96-well, rigid frame PC/PP, green, wells transparent 781544






L	127.76 mm
W	85.48 mm
H	16.1 mm
H _w	15.1 mm
X1	14.38 mm
Xn	113.38 mm
Y1	11.24 mm
Yn	74.24 mm
P	9.00 mm
D	5.46 mm
C	12-well
R	8-well

Well Number	96
Type	full skirt, low profile
Nominal Volume	150 µl
Material	PC/PP
Cut corner	H1
Quality	BIO-CERT® PCR QUALITY



Subject to technical modification without notice. Errors excepted.

PCR-Platte 96-well, rigid frame PC/PP, blue, wells transparent 781545







L	127.76 mm
W	85.48 mm
H	16.1 mm
H _w	15.1 mm
X1	14.38 mm
X _n	113.38 mm
Y1	11.24 mm
Y _n	74.24 mm
P	9.00 mm
D	5.46 mm
C	12-well
R	8-well

Well Number	96
Type	full skirt, low profile
Nominal Volume	150 µl
Material	PC/PP
Cut corner	H1
Quality	BIO-CERT® PCR QUALITY



Subject to technical modification without notice. Errors excepted.

PCR-Platte 96-well, rigid frame PC/PP, white, wells white 781546



L	127.76 mm
W	85.48 mm
H	16.1 mm
H _w	15.1 mm
X1	14.38 mm
X _n	113.38 mm
Y1	11.24 mm
Y _n	74.24 mm
P	9.00 mm
D	5.46 mm
C	12-well
R	8-well


Well Number	96
Type	full skirt, low profile
Nominal Volume	150 µl
Material	PC/PP
Cut corner	H1
Quality	BIO-CERT® PCR QUALITY



Subject to technical modification without notice. Errors excepted.

PCR-Platte 96-well, rigid frame, PC/PP, transparent wells transparent 781547

PRODUCT
INFORMATION




Well Number	96
Type	semi skirted, standard profile
Nominal Volume	200 µl
Material	PC/PP
Cut corner	A12
Quality	BIO-CERT® PCR QUALITY



Subject to technical modification without notice. Errors excepted.

PCR-Platte 96-well, rigid frame PC/PP, blue, wells transparent 781548





Well Number	96
Type	semi skirted, standard profile
Nominal Volume	200 µl
Material	PC/PP
Cut corner	A12
Quality	BIO-CERT® PCR QUALITY



Subject to technical modification without notice. Errors excepted.

PCR-Platte 96-well, rigid frame PC/PP, transparent wells white, Roche® Lightcycler® 781549



L	127.76 mm
W	85.48 mm
H	15.6 mm
H _w	15.1 mm
X1	14.38 mm
Xn	113.38 mm
Y1	11.24 mm
Yn	74.24 mm
P	9.00 mm
D	5.46 mm
C	12-well
R	8-well

Well Number	96
Type	semi skirted, low profile
Nominal Volume	150 µl
Material	PC/PP
Cut corner	H12
Quality	BIO-CERT® PCR QUALITY



Subject to technical modification without notice. Errors excepted.