



PCR plate 24-well 781411



L	81.0 mm
W	27.0 mm
H	20.9 mm
Hr	18.5 mm
P	9.0 mm
C	3 well
R	8 well



Type	non-skirted, standard
Well rim	not elevated
Volume	200 µl / well
Material	PP
Color	transparent
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Free of PCR inhibitors



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

Subject to technical modification without notice. Errors excepted.

PCR plate 24-well 781412



L	81.0 mm
W	27.0 mm
H	20.9 mm
Hr	18.5 mm
P	9.0 mm
C	3 well
R	8 well

Type	non-skirted, standard
Well rim	not elevated
Volume	200 µl / well
Material	PP
Color	white
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Free of PCR inhibitors



* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

Subject to technical modification without notice. Errors excepted.