


# Deep-well plate 1.1 ml 701352



L	127.50 mm
W	85.50 mm
H	41.00 mm
Hs	3.00 mm
X1	14.25 mm
Xn	113.25 mm
Y1	11.25 mm
Yn	74.25 mm
P	9.00 mm
D	7.20 mm
Hw	37.40 mm
C	12 well
R	8 well

<b>Well Number</b>	96
<b>Type</b>	Standard
<b>Bottom Shape</b>	U-bottom
<b>Nominal Volume</b>	1.1 ml
<b>Material</b>	PS
<b>Well Shape</b>	Round
<b>Cut corner</b>	H1
<b>Quality</b>	BIO-CERT® PCR QUALITY
<b>Suitable sealing mat</b>	701360



Subject to technical modification without notice. Errors excepted.