

# Colour up your Lab

SAY GOODBYE TO MONOTONY



**VITLAB**   
Competence in Labware

# Competence in Plastic Labware

## Practical helpers for diverse laboratory uses

Tired of dreary grey on grey?

Say goodbye to monotony! With its new product line VITLAB provides variety and sets colourful accents in your lab.

The transport, bottling and decanting of granulates, powders, pastes and liquids often requires special containers. VITLAB is one of the leading manufacturers of high quality plastic laboratory products that have become essential in up to date laboratory practice, e.g. in the food industry.

There are good reasons for this: superior break resistance means that the danger of injury is greatly reduced and that the products can be used for a long period of time. The utilization of carefully selected, high performance plastics ensures that VITLAB plastic products have an excellent chemical resistance and are light-weight for easy handling.

The high quality colours of the new product line contribute to a reliable distinction between different materials you use and support visual identification.

In addition, the contrast-rich colours are a highlight for any lab!



## Graduated beakers, coloured, PP



Made in Germany

Graduated beakers (size 500 ml) in four different colours. Stackable. Transparent, with easy readable, printed scale and stable handle. To preserve markings, cleaning at no higher than 60 °C is recommended. Suitable for food stuffs according to EC directive No. 10/2011.

Colour	Graduation ml	Height mm	Ø mm	PU	Cat. No.
blue	10	140	100	12	481942
yellow	10	140	100	12	481943
red	10	140	100	12	481944
green	10	140	100	12	481945
Set: blue, yellow, red, green (1 item each)	10	140	100	1	4811111



## Measuring scoops, coloured, PP



Measuring scoops (size: 100 ml) in seven different colours. Also suitable as weighing scoops. With precisely formed filling edge and comfortable, stable handle. Easily readable volume specification on the upper side of the handle. Suitable for food stuffs according to EC directive No. 10/2011.

Colour	PU	Cat. No.
white	12	39694
red	12	396940
grey	12	396943
black	12	396944
yellow	12	396946
blue	12	396950
green	12	396952
bright blue	12	396955
Set: white, red, grey, black, yellow, blue, green, bright blue (1 item each)	1	3961111



## Wash-bottles, coloured, PE-LD/PP

Narrow-mouth bottles, made from PE-LD. Available in four different colours to facilitate ready identification. Screw cap and spray tube made from PP. Precise spray jet and optimized medium backflow through the continuous and finely drawn tip of the spray tube.

Colour	Volume ml	Thread GL	Height* mm	Ø mm	PU	Cat. No.
red	250	25	135	58	5	132603
red	500	25	180	74	5	132703
green	250	25	135	58	5	132605
green	500	25	180	74	5	132705
yellow	250	25	135	58	5	132606
yellow	500	25	180	74	5	132706
blue	250	25	135	58	5	132608
blue	500	25	180	74	5	132708
Set: red, green, yellow, blue (1 item each)	500	25	180	74	1	1327111

\* Height without spray tube



# More information about VITLAB®?

---

We will be happy to send you  
more information about

- VITLAB® laboratory products for volume measurement
- VITLAB® laboratory products for saving and storing
- VITLAB® catalogue with complete product range

or

- the options for customized imprints on  
VITLAB® laboratory products.

Distributor:

**VITLAB GmbH** • Linus-Pauling-Str. 1 • 63762 Grossostheim • Germany  
tel: +49 (0) 6026 9 77 99-0 • fax: +49 (0) 6026 9 77 99-30 • e-mail: [info@vitlab.com](mailto:info@vitlab.com) • [www.vitlab.com](http://www.vitlab.com)