


Vacuum tubing transparent DN 15 mm

[Link to the website](#)

Product description

The transparent vacuum tubing offers maximum performance for demanding laboratory applications. Thanks to its transparency, it enables reliable visual control of the flow, while its excellent chemical resistance guarantees a wide range of applications. The smooth inner wall effectively prevents deposits and simplifies cleaning processes. The optimum balance of flexibility and stability enables both easy installation and a secure fit at connections - even in complex set-ups. Particularly noteworthy is the extremely long service life thanks to the excellent abrasion and bending resistance. With its reinforced construction, the hose offers maximum operational reliability for vacuum applications.

Price per meter, please choose desired length via quantity selection

Performance features

- Price per meter, please choose desired length via quantity selection
- Excellent chemical resistance
- Smooth inner wall prevents deposits and makes cleaning easier
- Long service life due to high abrasion and bending resistance
- Maximum operational safety for vacuum applications

Ordering data



Product name	Nominal mains voltage mains frequency	Plug type	Articleno.
DN 15 mm			20686063

Items supplied

Vacuum tubing transparent DN 15 mm (sale per meter).

Technical data

Dimension B in mm	34.9mm
Dimension C in mm	15.9mm

Dimension sheet

VACUUBRAND GMBH + CO KG
P.O. Box 1664
97877 Wertheim
T +49 9342 808 5550
F +49 9342 808 5555

info@vacuubrand.com
www.vacuubrand.com

VACUUBRAND®

VACUUBRAND is part of the Brand group, a manufacturer of high-quality and innovative laboratory equipment and vacuum technology with a global team of approx. 1,000 employees.

With highly motivated and qualified employees, we manufacture the world's most comprehensive product family for vacuum generation, measurement and control for the rough and fine vacuum range in the laboratory.

All with one goal:

Enabling our customers to run reliable, safe and efficient processes. This is why we commit ourselves to continuous innovation, first-class performance, and outstanding quality.

Engineering excellence is central to our ability to deliver products and solutions renowned for their quality, durability, and reliability. Innovative technology, high quality materials and extensive testing deliver unmatched performance - engineered and manufactured in Germany.

Together with you, our experts will find the best solution for your application and provide you with first-class service along the way.

© VACUUBRAND GMBH + CO KG

BRAND (Shanghai) Trading Co., Ltd.
Shanghai, China

T +86 21 6422 2318
info@brand.com.cn
china.brand.com.cn

BRAND Scientific Equipment Pvt. Ltd.
Mumbai, India

T +91 22 42957790
customersupport@brand.co.in
www.brand.co.in

BRANDTECH® Scientific, Inc.
Essex, CT. United States of America

T +1 860 767 2562
info@brandtech.com
www.brandtech.com