


Tuyau à vide, transparent DN 6 mm

[Lien vers le site web](#)

Description du produit

Le tuyau à vide transparent offre des performances maximales pour les applications de laboratoire exigeantes. Il permet un contrôle visuel fiable du flux, tandis que son excellente résistance chimique garantit des possibilités d'utilisation polyvalentes. La paroi intérieure lisse réduit efficacement les dépôts et simplifie le nettoyage. L'équilibre optimal entre flexibilité et stabilité permet à la fois une installation facile et une fixation sûre aux raccords, même dans les montages complexes. Il convient de souligner la durée de vie extrêmement longue grâce à l'excellente résistance à l'abrasion et à la flexion. Grâce à sa construction renforcée, le tuyau offre une sécurité de fonctionnement maximale pour les applications de vide. Prix par mètre, veuillez sélectionner la longueur souhaitée dans le choix de quantité

Caractéristiques

- Prix par mètre, veuillez sélectionner la longueur souhaitée dans le choix de quantité
- Excellente résistance chimique
- La paroi intérieure lisse empêche les dépôts et facilite le nettoyage
- Longue durée de vie grâce à une grande résistance à l'abrasion et à la flexion
- Sécurité de fonctionnement maximale pour les applications sous vide

Référence de commande



Désignation	Tension fréquence nominales	Type de fiche secteur	Référence
DN 6 mm			20686060

Livraison

Tuyau à vide transparent DN 6 mm (vente au mètre).

Caractéristiques techniques

Dimensions B in mm	15.9mm
Dimensions C in mm	6.4mm

Plan de masse

VACUUBRAND GMBH + CO KG
P.O. Box 1664
97877 Wertheim
T +49 9342 808 5550
F +49 9342 808 5555

info@vacuubrand.com
www.vacuubrand.com

VACUUBRAND®

VACUUBRAND is part of the Brand group, a manufacturer of high-quality and innovative laboratory equipment and vacuum technology with a global team of approx. 1,000 employees.

With highly motivated and qualified employees, we manufacture the world's most comprehensive product family for vacuum generation, measurement and control for the rough and fine vacuum range in the laboratory.

All with one goal:

Enabling our customers to run reliable, safe and efficient processes. This is why we commit ourselves to continuous innovation, first-class performance, and outstanding quality.

Engineering excellence is central to our ability to deliver products and solutions renowned for their quality, durability, and reliability. Innovative technology, high quality materials and extensive testing deliver unmatched performance - engineered and manufactured in Germany.

Together with you, our experts will find the best solution for your application and provide you with first-class service along the way.

© VACUUBRAND GMBH + CO KG

BRAND (Shanghai) Trading Co., Ltd.
Shanghai, China

T +86 21 6422 2318
info@brand.com.cn
china.brand.com.cn

BRAND Scientific Equipment Pvt. Ltd.
Mumbai, India

T +91 22 42957790
customersupport@brand.co.in
www.brand.co.in

BRANDTECH® Scientific, Inc.
Essex, CT. United States of America

T +1 860 767 2562
info@brandtech.com
www.brandtech.com